



Technical Note:

Preparing a Suspension Culture of Sf9 and Sf21 Insect Cells in a Spinner Flask

Introduction

A basic protocol is provided below to enable the transfer of adherent Sf9 and/or Sf21 cells from tissue culture flasks into a suspension culture in Bellco spinner flasks. This procedure allows for the initiation and maintenance of suspension cultures of Sf9 and Sf21 cells.

Materials required

- Adherent Sf9 or Sf21 insect cells
- Culture medium, with serum or serum free.
- Spinner flasks Bellco #1965 Series
- Magnetic stir plate
- 26–28°C incubator, CO₂ not required

Procedure

When initiating suspension cultures of Sf9 or Sf21 cells, log phase cells of 95% viability or greater should be used to ensure a successful initiation process.

To rapidly initiate a suspension culture of Sf9 or Sf21 cells:

- Grow up enough log phase adherent cells to start a spinner flask of the desired size with 0.5×10^6 viable cells/mL.
 - **Note:** Seeding at lower cell densities can cause increased doubling times during the initial stages of culture.
 - **Note:** Bellco recommends starting with a 100-mL or 250-mL spinner flask as these require fewer total cells than larger spinner flasks to initiate suspension cultures.
- Prepare clean, sterile spinner flask(s) with one of the following media types:
 - Grace's Insect Medium, Unsupplemented, plus 10% FBS.
 - Grace's Insect Medium, Supplemented, plus 10% FBS.
 - Serum Free Media
 - Gibco Sf-900 II SFM
 - Gibco Sf-900 III SFM
 - The recommended working volume of media is dependent on the size of the spinner flask, see below for recommended volumes based on maximizing aeration/oxygenation of the culture.
 - 100mL flask. Minimum working volume 30mL. Maximum working volume 50mL.
 - 250mL flask. Minimum working volume 80mL. Maximum working volume 125mL
 - 500mL flask. Minimum working volume 200mL. Maximum working volume 250mL
- Detach cells from the tissue culture flask(s) and count cells to determine cell density and viability. Cell viability needs to be >95% to ensure successful adaptation.

- **Note:** Cells can typically be easily detached using the sloughing method which dislodges the cell monolayer in the T-flask with the least manipulation and mechanical force which results in higher cell viabilities than other methods.
- **Basic sloughing method.**
 - Pour off to waste all but 5mL of media from the flask, irrespective of flask size.
 - Tilt the flask on its end so that all remaining media flows to one corner of the flask and away from the cells.
 - Draw some of the remaining media into a sterile pipette and stream the media across the cells starting at the bottom corner of the flask
 - It is recommended to use a side-to-side streaming motion and to move from the bottom corner up to the top opposite corner of the T-flask.
 - Always use a gentle media stream to dislodge the cells.
- Seed a clean, sterile spinner flask of the desired size with enough cells to reach a density of 1×10^6 cells/mL.
 - **Example:** A 100-mL spinner with 50 mL of culture at 0.5×10^6 viable cells/mL would require a total of $(0.5 \times 10^6 \text{ cells/mL}) \times (50 \text{ mL}) = 2.5 \times 10^7$ viable cells.
- Incubate the spinner flask(s) at 26–28°C with constant stirring at 80–90 rpm.
- When the cells reach a density of about 2×10^6 to 4×10^6 cells/mL, passage them at a seeding density of between 3×10^5 to 5×10^5 viable cells/mL into appropriately sized flask(s).

References

1. Bellco Technical Report TN-2025-001. General Recommendations for Using Bellco Spinner Flasks for Insect Cell Suspension Cultures.
2. Bellco Technical Report TN-2025-002. Key Points to Consider in Spinner Flask Cell Culture.
3. Gibco User guide MAN0000030. Growth and Maintenance of Insect Cell Lines. Rev C.0

Bellco Biotechnology

340 Edrudo Drive
 Vineland, NJ, 27006, USA
 Tel: +1 856 691 1075
 Fax: +1 856 691 3247
 Web: www.bellcoglass.com