

## **HELPFUL TIPS FOR HETTICH'S CENTRIFUGES:**

**The following questions could help you narrow down the correct centrifuge for your application;**

**~ Tube Type – Application?**

**Microliter Tubes**

**Round Bottom Tubes**

**Conical /Shirted Tubes**

**Blood Bags**

**High Capacity Tubes**

**Plates**

**~ Tube, Bag or Plate Size?**

**~ How many per run?**

**~ At what RCF?**

**~ Do you need refrigeration?**

**~ Do you need bio-containment?**

**~ Do you need a swing out or fixed rotor?**

**~ Cost / Budget constraints**

**Note: swing out rotors during centrifugation tubes are swung horizontal**

**Fixed angle rotors – used primarily in higher speed range**

**~ RPM – revolutions per minute**

