


**Magnetic Stirring Bars:**

Molded PTFE design with no seams

Excellent for mixing & buffer solutions

Ideal for procedures using all Bellco magnetic stirrers

PTFE Magnetic Stirr Bars



Stock No.	Description	"A"	"B"	
1975-00100	PTFE Mag Stirring Bar	7.5	25.4	
1975-00158	PTFE Mag Stirring Bar	9.5	41.3	
1975-00200	PTFE Mag Stirring Bar	9.5	50.8	
1975-00300	PTFE Mag Stirring Bar	12.1	76.2	
1975-00314	PTFE Mag Stirring Bar	19	82.6	
Dimensions are in (mm)				
Material:	PTFE Coating			

FYI - Information

F/G #	Blank Stock#	Usage
1975-00100	1975-0100A	1967-30025
1975-00158	1975-0158A	1965-30100,30500, 60001,60005,60015,80101,80251,80511,1969-30050
1975-00200	1975-0200A	1975-00200, & specials
1975-00300	1975-0300A	1964-03420,06420,1965-31000,60010,60014,60034,81001,83001 & 1969-31000
1975-00314	1975-0314A	1964-03410,06410,1965-31500,38000,60060,60150,60360 & 1969-31500

Magnet: Alnico #5