


SST Wire Culture Tube Racks:

All Bellco Culture Tube Racks are assembled with “wrap-around” welded wire construction to withstand rough handling. Racks are available for the following stock tube sizes; 13mm, 16mm, and 18mm, 20-22mm, 25mm & 38mm. They are manufactured from heavy gauge stainless steel wire. They will not rust or deteriorate in any way in the water bath or CO2 incubators. Packed one per case
Material: 303 SST wire with electropolish finish

Specifications of racks:

Tube diameter:	13-20mm	25mm & larger
Side Frame & uprights:	0.135”OD wire	0.162” OD wire
Cross Wire & Bottom mesh:	0.063” OD wire	0.063” OD wire



Tube Opening in Rack: “X”

For 13mm OD tubes:	9/16” (14.3mm)
For 16mm OD tubes:	11/16” (17.5mm)
For 18mm OD tubes:	13/16” (20.6mm)
For 20-22mm OD tubes:	15/16” (23.8mm)
For 25mm OD tubes:	1 1/16” (26.9mm)
For 38mm OD tubes:	1 9/16” (39.6mm)