



# ROLLER APPARATUS

Keep your research rolling



Manufactured in  
Vineland, NJ  
since 1936

**MADE IN** ★ ★ ★  
**USA** ★ ★ ★



## ROLL WITH PRECISION. GROW WITH CONFIDENCE.

Bellco's Roller Apparatus fits seamlessly in incubators and scales with your workflow. With customizable spacing, non-slip chemically resistant rollers, digital touchscreen control, and connectivity via Wi-Fi or Ethernet, it's precision-made for high-yield cell culture.

## DESIGNED FOR HOW YOU WORK

Roller Apparatus can be designed and engineered to fit your needs and requirements. Just ask us!

### MODELS

Modular (easy to expand)  
Fixed Frame

### NUMBER OF BOTTLES

1 to 176  
Custom Quantities Available

### ROLLER BOTTLE VOLUMES

670, 850, 1170, 1330, & 1585mL

### COMMUNICATION

Ethernet  
Wi-Fi  
4-20mA  
RS-485

### LOCATION

Bench Top  
Floor Mount  
with Casters

### VOLTAGE

100V  
115V  
230V

### CONTROLLER

Digital Display  
PID with Digital Display  
PID with Touchscreen Interface

### OPTIONS

Battery Back-up  
Temperature Display  
Motion Alarm  
Data Logger

### NUMBER OF DECKS

1 to 11 position

### DECK SPACING

6.375in  
7.125in  
Custom Spacing

### ROLLERS

Chemically Resistant  
Non-Slip EDPM

### RPM

0.1 to 4.0rpm  
3.0 to 9.0rpm  
Custom Settings

### DRIVES LOCATIONS


Top Mount  
Bottom Mount  
Side Mount





Made for each other.

## Bellco Incubators and Roller Apparatus.



Bellco's Roll-In Incubators are built to accommodate large-scale roller apparatus setups with ease. Designed for uniform temperature distribution and optimal airflow, these units maintain a consistent environment for high-yield cell culture. With durable stainless-steel interiors, easy-access doors, and reliable controls, Bellco incubators make scaling up as simple as rolling in.



SCAN ME  
To learn more!



**CONTACT US:**

**Phone:** 1-800-257-7043

**Email:** [cservice@bellcoglass.com](mailto:cservice@bellcoglass.com)

**Web:** [www.bellcoglass.com](http://www.bellcoglass.com)

