



# BRUSHLESS DIGITAL ROLLER DRUM DRIVE

## FLEXIBLE SPEED. HIGH-CAPACITY PERFORMANCE.

Supports a broad range of applications, from blood processing to viral cultures, with customizable speed settings and drum capacities to match your lab's workflow.

One size  
fits most.  
**CUSTOM**  
fits the rest.

**DRUM TILT**  
7°

**DIGITAL DISPLAY**

**EASY MOUNTING BRACKET**

### MODELS

Standard Brush-less, Low Speed  
Standard Brush-less  
Low Profile Brush-less, Low Speed  
Low Profile Brush-less

**AVAILABLE  
VOLTAGE**  
115V  
230V

**SPEED RANGE**  
2 - 60 rpm  
0.5 - 15 rpm  
Custom Speeds  
Available



**MADE IN** ★★  
**USA** ★★

\*Most Bellco products can be configured to suit your unique needs.

# ROLLER DRUMS



Part # 7736-20025

**156 POSITION  
25mm TUBES**



Part # 7736-20050

**128 POSITION  
30mm TUBES**



Part # 7736-20164

**164 POSITION  
16mm TUBES**



Part # 7736-20351

**351 POSITION  
16mm TUBES**



Part # 7736-21106

**106 POSITION  
LOW PROFILE  
18mm TUBES**

## CONTACT US:

**Phone:** 1-800-257-7043

**Email:** [cservice@bellcoglass.com](mailto:cservice@bellcoglass.com)

**Web:** [www.bellcoglass.com](http://www.bellcoglass.com)