



# DIGITAL ORBITAL SHAKERS

## COMPACT POWER. CUSTOM FLEXIBILITY.

Versatile speed, swappable trays, and timed control,  
all in one durable, space-saving shaker.

One size  
fits most.  
CUSTOM  
fits the rest.

### OPTIONS

Standard  
Mini

### AVAILABLE VOLTAGE

100V  
115V  
230V

### TRAY OPTIONS

10" x 10" (25 x 25 cm) Tray  
20" X 20" (51 x 51 cm) Tray  
Mini-Shaker Tray 8-Position  
Mini-Shaker Tray 32-Position

### MAXIMUM LOAD

24 lbs (11 kg)

### SPEED RANGE

0 – 350 rpm  
40 – 1550 rpm

### ORBIT

1" (25 mm)  
1/16" (1.6 mm) - Mini



**MADE IN** ★ ★ ★  
**USA** ★ ★ ★

\*Most Bellco products can be configured to suit your unique needs.



# DIGITAL ORBITAL SHAKERS

## DIGITAL ORBITAL SHAKER

The Bellco Digital Shaker is a compact, heavy-duty unit for agitating gels and assay media, perfect for staining and destaining. It offers variable speeds from 30–350 rpm and supports a range of vessels with interchangeable trays, including flask clamp and peg-and-grommet styles. A timer allows automatic shutoff from 2–120 minutes.



### Specifications:

- Speed Range: 30 – 350 rpm +/- 2%
- Circular Orbit: 1" (25 mm)
- Base Overall Dimensions:  
12" x 14" x 4 1/2"  
(13 x 28 x 28 cm)

## DIGITAL MINI ORBITAL SHAKER

The Bellco Digital Mini-Orbital Shaker is a high speed shaker for general laboratory cell culture experimentation. The design incorporates an orbital platform that accepts a variety of optional trays. Peg and grommet design offers a simple method to install or change 4 versions of trays. It has a brushless DC motor designed for minimal maintenance.



### Specifications:

- Speed Range: 40 – 1550 rpm
- Circular Orbit: 1/16" (1.6 mm)
- Base Overall Dimensions:  
13 3/4" x 14 1/2" x 5 1/2"  
(35 x 37 x 14 cm)

### CONTACT US:

**Phone:** 1-800-257-7043

**Email:** [cservice@bellcoglass.com](mailto:cservice@bellcoglass.com)

**Web:** [www.bellcoglass.com](http://www.bellcoglass.com)

