



HIGH-TORQUE OVERHEAD DRIVES

HIGH TORQUE. LOW MAINTENANCE.

Engineered for high-viscosity and large-volume mixing, Bellco's Digital Overhead Drive ensures controlled, consistent, and aseptic performance.

One size
fits most.
CUSTOM
fits the rest.

DISPLAY

Digital Display
Analog
Direct Drive

REVERSE DIRECTION SWITCH

AVAILABLE VOLTAGE

115V
230V

TORQUE

Hi Torque 31 In-Lb Torque

MODELS

High Torque
Digital Control
High Torque
Process Control

SPEED RANGE

5 – 150 rpm

EASY MOUNTING BRACKET
WITH 'PEG & GROMMET' DESIGN



MADE IN ★★
USA ★★

*Most Bellco products can be configured to suit your unique needs.



HIGH-TORQUE OVERHEAD DRIVES

HIGH-TORQUE DIGITAL OVERHEAD DRIVE

Engineered for cell culture and aseptic mixing, the Bellco High Torque Digital Overhead Drive delivers powerful, direct-drive agitation, perfect for medium to high-viscosity liquids. With a quiet, brushless DC motor and CE-compliant 24VDC power, it offers reliable, maintenance-free performance for demanding lab environments.



Specifications:

- Speed Range: 5 – 150 Rpm +/- 2%
- Brushless Motor For Low Maintenance
- Reverse Direction Switch
- Dimensions: 9.5" X 6" X 7.5"
(24 X 15 X 19cm) Includes Handle
- Hi-Torque Models: 31 In-Lb Torque, 50w

HIGH-TORQUE PROCESS CONTROLLED OVERHEAD DRIVE

Designed for cell culture and aseptic mixing, the Bellco High Torque Overhead Drive delivers powerful, direct-drive agitation without magnetic decoupling, ideal for high-viscosity or large-volume applications. Featuring a quiet, brushless DC motor and CE-compliant 24V power, it's a reliable, maintenance-free solution for demanding



Specifications:

- Speed Range: 5 – 150 Rpm +/- 1%
- Brushless Motor For Low Maintenance
- Reverse Direction Switch
- Dimensions: 9.5" X 6.0" X 7.5"
(24 X 15 X 19cm) Includes Handle
- Hi-Torque Models: 31 In-Lb Torque, 50w

CONTACT US:

Phone: 1-800-257-7043

Email: cservice@bellcoglass.com

Web: www.bellcoglass.com

