



BENCHTOP ORBITAL SHAKERS

COMPACT POWER. CUSTOM FLEXIBILITY.

Versatile speed, interchangeable trays, and timed control,
all in one durable, space-saving shaker.

One size
fits most.
CUSTOM
fits the rest.

TRAY OPTIONS

10 x 10 inch (25 x 25 cm) Tray
10 x 14 inch (25 x 35 cm) Tray
12 x 12 inch (30 x 30 cm) Tray

20 x 20 inch (51 x 51 cm) Tray
Mini-Shaker Tray 8-Position
Mini-Shaker Tray 32-Position

OPTIONS

Standard
Mini

**MAXIMUM
LOAD**
24 lbs

AVAILABLE VOLTAGE

100V
115V
230V

ORBIT

1"
1/16" (Mini)

SPEED RANGE

30 – 350 rpm
50 – 1550 rpm (Mini)



MADE IN ★★
USA ★★

*Most Bellco products can be configured to suit your unique needs.



BENCHTOP ORBITAL SHAKERS

BENCHTOP ORBITAL SHAKER

Bellco Bench-top Orbital Shakers offer reliable mixing for gel staining, washing, and general agitation. With smooth speed control from 30 to 350 rpm ($\pm 5\%$), they deliver consistent performance across a wide range of lab applications. Available in two platform sizes: 25 x 25 cm and 51 x 51 cm, to suit your workflow needs.



Specifications:

- Speed Range: 30 – 350 rpm $\pm 5\%$
- Circular Orbit: 1" (2.5 cm)
- Maximum Load: 11 Kg (25 Lbs)
- Dimensions: 5" x 11" x 11"
(13.0 x 28 x 28 cm)

BENCHTOP MINI ORBITAL SHAKER

Bellco Mini Orbital Shakers deliver consistent, true orbital motion for optimal mixing of 96-well micro tube plates in cell culture, virology, immunology, and serology. With a broad speed range of 50–1550 rpm, timed or continuous operation, and a cast iron base to prevent walking, these compact shakers are built for stability and reproducibility.



Specifications:

- Speed Range: 50 – 1550 rpm $\pm 5\%$
- Circular Orbit: 1/16" (1.6 mm)
- Base Overall Dimensions:
5" x 11" x 11"
(13.0 x 28 x 28 cm)

CONTACT US:

Phone: 1-800-257-7043

Email: cservice@bellcoglass.com

Web: www.bellcoglass.com

