



OVERHEAD DRIVES

BUILT FOR ROUTINE MIXING. DESIGNED TO LAST.

Engineered for high-viscosity and large-volume mixing, Bellco's Digital Overhead Drive ensures controlled, consistent, and aseptic performance.

One size
fits most.
CUSTOM
fits the rest.

SPEED RANGE

5 – 150 rpm

REVERSE DIRECTION SWITCH

AVAILABLE VOLTAGE

115V
230V

TORQUE

Standard 19 In-Lb Torque

MODELS

Digital Control
Process Control
High Torque
Digital Control
High Torque
Process Control

DISPLAY

Digital Display
Analog
Direct Drive



EASY MOUNTING BRACKET
WITH 'PEG & GROMMET' DESIGN

MADE IN ★ ★ ★
USA ★ ★ ★

*Most Bellco products can be configured to suit your unique needs.

DIGITAL OVERHEAD DRIVE

Engineered for reliable, direct-drive mixing in viscous or large-volume applications. A quiet, brushless DC motor ensures maintenance-free operation, while CE compliance and a 24VDC power supply support aseptic performance in demanding lab environments.



Specifications:

- Speed Range: 5 – 150 rpm +/- 2%
- Brushless Motor For Low Maintenance
- Reverse Direction Switch
- Easy Mounting Bracket With 'Peg & Grommet' Design
- Dimensions: 4.5" X 5.5" X 8" (11 X 14 X 20 cm) Includes Handle

PROCESS CONTROLLED OVERHEAD DRIVE

Engineered for reliable, high-volume mixing, the Bellco Overhead Drive eliminates magnetic decoupling issues in viscous or large-volume applications. Featuring a quiet, brushless DC motor and maintenance-free design, it delivers powerful, aseptic agitation with CE-compliant performance.



Specifications:

- Speed Range: 5 – 150 rpm +/- 1%
- Brushless Motor For Low Maintenance
- Reverse Direction Switch
- Easy Mounting Bracket With 'Peg & Grommet' Design
- Dimensions: 4.5" X 5.5" X 8" (20 X 12 X 17 cm) Includes Handle
- Standard Models: 19 In-Lb Torque, 30w

CONTACT US:

Phone: 1-800-257-7043

Email: cservice@bellcoglass.com

Web: www.bellcoglass.com