



DIGITAL ROCKER

CONSISTENT MOTION, RELIABLE RESULTS

Engineered for gentle, uniform rocking, Bellco platforms support optimal cell growth and adaptable workflows.

One size
fits most.
CUSTOM
fits the rest.

TRAY OPTIONS

25 x 25 cm Tray
51 x 51 cm Tray
3-Tier Rocker Tray

**MAXIMUM
LOAD**
9 kg

**RANGE OF
MOTION**
0 - 12 Degrees
(Adjustable)

SPEED RANGE
0 - 30
Oscillations
per Minute
±2%

CONTROLLER
Digital

TIMER RANGE
3 - 120 Minutes



MADE IN ★★
USA ★★

*Most Bellco products can be configured to suit your unique needs.



STANDARD ROCKER

CONSISTENT MOTION, RELIABLE RESULTS

Engineered for gentle, uniform rocking, Bellco platforms support optimal cell growth and adaptable workflows.

TRAY OPTIONS

25" x 25" Tray
51" x 51" Tray
Three Tier Rocker Tray
Mini-Shaker Tray 8-Position
Mini-Shaker Tray 32-Position

MAXIMUM

LOAD

20 lbs
(9 kg)

RANGE OF MOTION

7 Degrees
(Fixed)

SPEED RANGE

0 - 12
Oscillations
per Minute



CONTROLLER

Analog

VOLTAGE

100V, 115V, OR 230V

CONTACT US:

Phone: 1-800-257-7043

Email: cservice@bellcoglass.com

Web: www.bellcoglass.com

