



2D BARCODED GLASS VIALS

Bellco Glass introduces our 9mm Borosilicate Glass 2D Barcoded Vials, engineered for durability and precision. Made to withstand harsh laboratory solvents, abrasive conditions, and a wide temperature range from -80°C to 150°C, these vials are ideal for demanding research environments.

Key Features:

Chemical & Thermal Resistance

- Permanent labels that are abrasion and solvent resistant.

Pre-Marked Convenience

- Ready-to-use vials eliminate manual labeling errors.

Contamination Control

- Borosilicate glass lowers leaching risk along with least pH shift.

Seamless Integration

- Optimized for automation and inventory systems.

Various Capacities

- Vials available in 1mL, 1.5mL, 2mL

Automation Friendly

- The Multi-Tier Base Plate follows ANSI footprint

Applications

- Pharmaceuticals
- Environmental Testing
- Biotech & Research
- Chromatography (HPLC & GC)
- Forensics





2D BARCODED GLASS VIALS

Description	Case Qty
1.0ml vial 9x30mm bottom barcoded loaded in base plate	5 racks of 96
1.0ml vial 9x30mm bottom & side barcoded loaded in base plate	5 racks of 96
1.0ml vial 9x30mm bottom barcoded in vial loader	5 Vial Loaders
1.0ml vial 9x30mm bottom & side barcoded in vial loader	5 Vial Loaders
1.0ml vial 9x30mm bottom barcoded	100 vials
1.0ml vial 9x30mm bottom & side barcoded	100 vials
1.5ml vial 9mmx44mm bottom barcoded loaded in base plate	5 racks of 96
1.5ml vial 9mmx44mm bottom & side barcoded loaded in base plate	5 racks of 96
1.5ml vial 9x44mm bottom barcoded in vial loader	5 Vial Loaders
1.5ml vial 9x50mm bottom & side barcoded in vial loader	5 Vial Loaders
1.5ml vial 9x44mm bottom barcoded	100 vials
1.5ml vial 9x44mm bottom & side barcoded	100 vials
2.0ml vial 9mmx50mm bottom barcoded loaded in base plate	5 racks of 96
2.0ml vial 9mmx50mm bottom & side barcoded loaded in base plate	5 racks of 96
2.0ml vial 9x50mm bottom barcoded in vial loader	5 Vial Loaders
2.0ml vial 9x50mm bottom & side barcoded in vial loader	5 Vial Loaders
2.0ml vial 9x50mm bottom barcoded	100 vials
2.0ml vials 9x50mm bottom & side barcoded	100 vials



96 Round Well Clear Sealing Mat with Premium Silicone	5
96 Round Well Royal Blue Sealing Mat Premium Silicone	5
96 Round Well Pre-Slit Clear Sealing Mat with Premium Silicone	5
96 Round Pre-Slit Royal Blue Sealing Mat with Premium Silicone	5
Molded Tan PTFE/Silicone Plug - Single Plug	100
Molded Tan PTFE/Silicone Strip - 8 Plugs	12
Molded Tan PTFE/Silicone Strip - 12 Plugs	8
96 Round 9mm Blue TPE Cap Mat	10

CONTACT US:

Phone: 1-800-257-7043

Email: cservice@bellcoglass.com

Web: www.bellcoglass.com

