

**Stoppers:**

Guide for selecting rubber stoppers; (dimensions are in mm)

Stopper #	Top Diameter	Bottom Diameter	Length	For OD Tubes/Flasks	or	For ID Tubes/Flasks
000	12.7	8.2	21	3		----
00	15	10	25	12-15		10-13
0	17	13	25	16-18		13-15
1	19	14	25	19-20		15-17
2	20	16	25	20-21		16-18.5
3	24	18	25	22-24		18-21
4	26	20	25	25-26		20-23
5	27	23	25	27-28		23-25
5 ½	28	24	25	28-29		25-26
6	32	26	25	29-30		26-27
6 ½	34	27	25	30-34		27-31.5
7	37	30	25	35-38		30-34
8	41	33	25	38-41		33-37
9	45	37	25	41-45		37-41
10	50	42	25	45-50		42-46
10 ½	53	45	25	48-51		45-47
11	56	48	25	52-56		48-51.5
11 ½	63	50	25	57-61		51-56
12	64	54	25	62-64		54-59
13	68	58	25	64-68		58-53
14	90	75	39	80-90		75-85
15	103	83	39	92-100		83-95