



CARBOYS, ASPIRATOR & RESERVOIR BOTTLES

DESIGNED FOR HOW YOU WORK

Carboys, Aspirator & Reservoir Bottles can be tailored and engineered to fit your needs and requirements.

NECK AND CLOSURES

Stopper Neck, #12
Other sizes available
Threaded Neck
45mm, 70mm, or 100mm with Cap

VESSEL OPTIONS

Clear Glass
Amber Coated
Plastic Safety Coating

PACKAGING OPTIONS

Clean Packed,
Non-Sterile
Sterile
Depyrogenated

SIZE OPTIONS

1L to 150L

GRADUATION OPTIONS

None
Liters or Gallons
Select Increments

CONNECTION OPTIONS

Hose Connections
Side Arms

POSITION OPTIONS

Neck
Side Wall
Bottom Mounted

SCREENING OR SERIALIZATION

Available on Request





Precision-Engineered Customization for Superior Research Outcomes.

Our carboy aspirators and reservoir bottles deliver extensive customization options that directly enhance your laboratory's efficiency and research accuracy. When you can select the exact specifications your protocols demand, you eliminate compatibility issues, reduce contamination risks, and streamline your workflows.

Material Selection and Right-Sized Solutions

All Bellco Carboys, Aspirator and Reservoir bottle are manufactured from Type 1 Class A borosilicate glass and can be provided either as clear glass or amber glass for light sensitive products and can be plastic coated for added safety.

Right-Sized Solutions Prevent Waste and Optimize Storage.

With sizes from 1 to 150 liters, you can match vessel capacity precisely to your batch requirements - eliminating excess headspace that can compromise sterility, reducing material waste, and maximizing valuable lab space.



CONTACT US:

Phone: 1-856-691-1075

Email: cservice@bellcoglass.com

Web: www.bellcoglass.com

