



CARBOYS, ASPIRATOR & RESERVOIR BOTTLES

Type I Class A Borosilicate Glass
Made in the USA Since 1936



SIZE RANGE

1L to 150L Capacities
Standard & Custom Configurations



CLOSURE OPTIONS

#12 Stopper
Threaded Neck
Screw Cap



CONNECTION OPTIONS

Hose Connections
Side Arms
Bottom Outlets



GRADUATION OPTIONS

None
Liters
Gallons
Custom Increments



GLASS MODIFICATIONS

Clear
Amber Stained
BellCoat Safety Coated



TRACEABILITY & PACKAGING

2D Barcoding
Serialization
Clean Packed
Sterile

ONE SIZE FITS MOST. CUSTOM FITS THE REST.

From ports and hose connections to safety coatings and custom configurations, Bellco vessels can be tailored to your process without the cost or delays of new tooling.

PRODUCT SELECTION GUIDE


Standard Sizes & Configurations



CARBOYS

SKU	Volume	Neck Size	Graduated
5057-00005	9.5L	#12	
5057-00007	9.5L	#12	
5057-00002	13L	#12	
5057-00003	19L	#12	
5057-00009	19L	100mm	3L-15L
5057-00001	45L	#12	
5057-00011	45L	#12	5L-40L

COMPATIBLE CLOSURES



SKU	Description
2048-12000	#12 Silicone Stopper
A12-418	70mm PP Screw Cap, Blue
A523-146	100mm PP Screw Cap, Blue



RESERVOIR BOTTLES

SKU	Volume	Neck Size	Connection	Outlet	# Outlets
5057-00010	45L	#12	Hose Connection	5/8" hose barb	2



ASPIRATOR BOTTLES

SKU	Volume	Neck Size	Connection	Outlet	# Outlets	BellCoat
5057-00004	20L	#12	Hose Connection	3/8" hose barb	2	-
5057-00008	20L	#12	Hose Connection	15mm	1	-
5057-00101	1L	#12	Hose Connection	3/8" hose barb	1	Yes
5057-00102	2L	#12	Hose Connection	3/8" hose barb	1	Yes
5057-00106	5L	#12	Hose Connection	3/8" hose barb	1	Yes
5057-00107	10L	#12	Hose Connection	3/8" hose barb	1	Yes
5057-00109	20L	#12	Hose Connection	3/8" hose barb	1	Yes

GLASSWARE MODIFICATIONS



BellCoat Safety Coating

Helps contain glass fragments and contents if breakage occurs while improving grip & handling.



Amber Staining

Protects light- and UV-sensitive materials during storage, processing & transport.



2D Barcoding

Supports product identification, traceability & inventory management.



Bellco Glass, Inc.

1-800-257-7043 | cservice@bellcoglass.com | bellcoglass.com



Made in the USA Since 1936

Scan for Product
Specifications &
Ordering Information

